

ADAPTERS




Adapters, cone/rubber tubing, right angled connection.

Code No.	Cone size	O/D nipple mm	Std. pack qty.
ATA/0	10/19	9	2
ATA/1	14/23	9	10
ATA/2	19/26	9	10
ATA/3	24/29	9	10
ATA/4	29/32	9	10
ATA/5	34/35	9	2




Adapters, cone/rubber tubing, straight connection.

Code No.	Cone size	O/D nipple mm	Std. pack qty.
ATS/1	14/23	9	10
ATS/2	19/26	9	10
ATS/3	24/29	9	10




Adapters, conical/spherical

Code No.	Cone size	Cup size	Std. pack qty.
ACC/2	19/26	S 29	5
ACC/3	24/29	S 35	5



Adapters, conical/spherical


Code No.	Socket size	Ball size	Std. pack qty.
ASB/2	19/26	S 29	2
ASB/3	24/29	S 35	2



Adapters, expansion.


Code No.	Socket size	Cone size	Std. pack Qty.
AX10	14/23	10/19	10
AX21	19/26	14/23	10
AX31	24/29	14/23	10
AX32	24/29	19/26	10
AX41	29/32	14/23	10
AX42	29/32	19/26	10
AX43	29/32	24/29	10
AX52	34/35	19/26	10
AX53	34/35	24/29	10
AX54	34/35	29/32	2
AX64	40/38	29/32	1

ADAPTERS




Adapters, multiple, two necks parallel.

Code No.	Socket size	Cone size	Std. pack qty.
AMP/00	10/19	10/13	2
AMP/11	14/23	14/23	5
AMP/22	19/26	19/26	5
AMP/32	19/26	24/29	5
AMP/53	19/26	34/35	2
AMP/33	24/29	24/29	2
AMP/44	29/32	29/32	2




Adapters, multiple, one parallel and one 45° neck.

Code No.	Socket size	Cone size	Std. pack qty.
AMA/11	14/23	14/23	2
AMA/22	19/26	19/26	2
AMA/32	19/26	24/29	2
AMA/33	24/29	24/29	2
AMA/43	24/29	29/32	2
AMA/53	24/29	34/35	2





Adapters, multiple, two parallel and one 45° neck.

Code No.	Sockets size Parallel	Angled	Cone size	Std. pack qty.
AM/321	19/26	14/23	24/29	2
AM/322	19/26	19/26	24/29	2
AM/422	19/26	19/26	29/32	1
AM/522	19/26	19/26	34/35	1




Adapter, receiver, intermediate (Perkin).

Code No.	Socket size	Cone size	Min. capacity ml	Std. pack qty.
RIP/11	14/23	14/23	10	1
RIP/22	19/26	19/26	32	1
RIP/23	19/26	24/29	32	1
RIP/33	24/29	24/29	32	1



Adapter, receiver, vacuum (to take 4 x 100 ml flasks).


Code No.	Socket size	Cone size	Std. pack qty.
ARA/21	19/26	14/23	2
Vertical delivery adapter ARA/22	19/26	19/26	2
Sloping delivery adapter ARA/22B	19/26	19/26	2



Adapter, receiver, vertical with vacuum connection.


Code No.	Socket size	Cone size	Std. pack qty.
RASV/33	24/29	24/29	2
RASV/22	19/26	19/26	2

ADAPTERS




Adapters, reduction.

Code No.	Socket size	Cone size	Std. pack qty.
AR01	10/19	14/23	5
AR12	14/23	19/26	20
AR13	14/23	24/29	20
AR14	14/23	29/32	20
AR23	19/26	24/29	20
AR24	19/26	29/32	5
AR25	19/26	34/35	5
AR26	19/26	40/38	5
AR27	19/26	45/40	5
AR29	19/26	55/44	5
AR34	24/29	29/32	5
AR35	24/29	34/35	5
AR36	24/29	40/38	5
AR37	24/29	45/40	5
AR38	24/29	50/42	5
AR39	24/29	55/44	5
AR45	29/32	34/35	5
AR47	29/32	45/40	5
AR56	34/35	40/38	5
AR57	34/35	45/40	5
AR58	34/35	50/42	5
AR59	34/35	55/44	5




Adapters, cone/screw thread.

Code No.	Cone size	Thread size	Diameter through Cap, mm	Std. pack qty.
AST/1	14/23	13	6.0-7.0	20
AST/2	19/26	13	6.0-7.0	20
AST/3	24/29	13	6.0-7.0	20
AST/4	29/32	13	6.0-7.0	20
AST/5	34/35	13	6.0-7.0	2
AST/6	14/23	18	7.0-8.5	20
AST/7	19/26	18	7.0-8.5	20
AST/8	24/29	18	7.0-8.5	20
AST/9	29/32	18	7.0-8.5	20
AST/10	34/35	18	7.0-8.5	2
AST/11	19/26	24	10.5-11.5	20
AST/12	24/29	24	10.5-11.5	20
AST/13	29/32	24	10.5-11.5	20
AST/14	34/35	24	10.5-11.5	2
AST/15	29/32	28	18.0-19.0	2
AST/16	34/35	28	18.0-19.0	2



Adapters, cone and two screwthreads. (one vertical and one at 90°).

Code No.	Cone size	Thread size	Diameter through Cap, mm	Std. pack qty.
ACS/1	14/23	13	6.0-7.0	1
ACS/2	14/23	18	7.0-8.5	1
ACS/3	19/26	18	7.0-8.5	1
ACS/4	24/29	24	10.5-11.5	1
ACS/5	29/32	24	10.5-11.5	1
ACS/6	29/32	28	18.0-19.0	1




ADAPTERS

Adapters, cone and three screwthreads.


Code No.	Cone size	Thread size vertical 90°		Std. pack qty.
		45°	18	
ACS/7	19/26	18	18	1
ACS/8	29/32	28	18	1
ACS/9	34/35	28	18	1

Adapters, socket/cone with 'T' connection.




Code No.	Socket size	Cone size	Std. pack qty.
ATR/11	14/23	14/23	5
ATR/12	14/23	19/26	5
ATR/13	14/23	24/29	5
ATR/22	19/26	19/26	5
ATR/23	19/26	24/29	5
ATR/24	19/26	29/32	2
ATR/25	19/26	34/35	2
ATR/33	24/29	24/29	5
ATR/44	29/32	29/32	2

Adapters, socket/rubber tubing, straight connection.



Code No.	Socket size	O/D nipple mm	Std. pack qty.
ATF/2	19/26	9	5
ATF/3	24/29	9	5
ATF/4	29/32	9	5

Adapters, stopcock, cone/rubber tubing, straight connection, interchangeable key.




Code No.	Cone size	O/D nipple mm	Key bore mm	Std. pack qty.
ASS/1	14/23	9	4	5
ASS/2	19/26	9	4	5
ASS/3	24/29	9	4	5
ASS/4	29/32	9	4	5

Adapters, stopcock, cone/rubber tubing, right angle connection, interchangeable key.




Code No.	Cone size	O/D nipple mm	Key bore mm	Std. pack qty.
ASA/1	14/23	9	4	5
ASA/2	19/26	9	4	5
ASA/3	24/29	9	4	5
ASA/5	34/35	9	4	5



ADAPTERS

Adapters, stopcock, socket/rubber tubing, interchangeable key.




<i>Code No.</i>	<i>Socket size</i>	<i>O/D nipple mm</i>	<i>Key bore mm</i>	<i>Std. pack qty.</i>
ASF1	14/23	9	4	5
ASF2	19/26	9	4	5
ASF3	24/29	9	4	5
ASF2B	19/26	9	4	2



Adapter, straight delivery.

<i>Code No.</i>	<i>Socket size</i>	<i>Overall length mm</i>	<i>Std. pack qty.</i>
RAS1	14/23	180	2
RAS2	19/26	180	5
RAS3	24/29	180	2
RAS4	29/32	180	2
RAS5	34/35	180	2



Adapters, swan neck (screwthread for use with thermometer or air leak tube).

<i>Code No.</i>	<i>Socket size</i>	<i>Cone size</i>	<i>Thread size</i>	<i>Std. pack qty.</i>
ASN1	14/23	14/23	13	2
ASN2	19/26	19/26	13	2
ASN3	19/26	24/29	13	2
ASN4	24/29	24/29	13	2