

Dean & Stark Apparatus.

Code No.	Capacity ml	Graduations ml	Socket size	Cone size	Std. pack qty.
Receiver					
DS1	2.0	0.05	19/26	24/29	2
DS2	7.5	0.10	19/26	24/29	2
DS3	25.0	0.20	19/26	24/29	2

Receiver with stopcock

DS4	10.00	0.10	19/26	24/29	2
-----	-------	------	-------	-------	---



Dean & Stark Apparatus, modified, heavy entrainer.

Code No.	Capacity ml	Graduations ml	Socket size	Cone size	Std. pack qty.
Receiver					
DS5	3.0	0.10	19/26	24/29	2
DS6	12.5	0.10	19/26	24/29	1



Dephlegmator.

Code No.	Cone size	Effect. length mm	Total dia. mm	Finger dia. mm	Std. pack qty.
CCF14/00	19/26	100	175	12	5

DISTILLATION

Distillation columns, plain.


Code No.	Socket size	Cone size	Effective length mm	Std. pack qty.
FCP 25/33	24/29	24/29	250	2
FCP 50/33	24/29	24/29	500	2
FCP 50/44	29/32	29/32	500	1
FCP 75/33	24/29	24/29	750	1
FCP 250/53	34/35	24/29	250	2
FCP 500/53	34/35	24/29	500	1

Distillation column, Vigreux.

Code No.	Socket size	Cone size	Effective length mm	Std. pack qty.
FCV 36/33	24/29	24/29	360	2
FCV 60/33	24/29	24/29	600	1

Distillation heads, Claisen, sloping.

Code No.	Cone to fit flask	Cone to fit condenser	Sockets size	Std. pack qty.
SBC/11	14/23	14/23	14/23	5
SBC/22	19/26	19/26	14/23	5
SBC/23	24/29	19/26	14/23	5
SBC/33	24/29	24/29	14/23	1
SBC/44	29/32	29/32	14/23	1
SBC/25	34/35	19/26	14/23	1




DISTILLATION HEADS




Distillation, splash head, adapter.

<i>Code No.</i>	<i>Socket size</i>	<i>Cone size</i>	<i>Std. pack qty.</i>
SPA/23	19/26	24/29	2
SPA/33	24/29	24/29	2




Distillation heads, steam distillation, sloping.

<i>Code No.</i>	<i>Cone to fit flask</i>	<i>Cone to fit condenser</i>	<i>Std. pack qty.</i>
SPD/23	24/29	19/26	5
SPD/25	34/35	19/26	1



Distillation heads, plain with thermometer socket.

<i>Code No.</i>	<i>Cone to fit flask</i>	<i>Cone to fit condenser</i>	<i>Socket size</i>	<i>Std. pack qty.</i>
SBS/00	10/13	10/13	10/19	5
SBS/11	14/23	14/23	14/23	10
SBS/12	19/26	19/26	14/23	10
SBS/22	19/26	19/26	19/26	5
SBS/23	24/29	19/26	14/23	5
SBS/33	24/29	24/29	14/23	5
SBS/24	29/32	19/26	14/23	1
SBS/44	29/32	29/32	14/23	2
SBS/25	34/35	19/26	14/23	1
SBS/35	34/35	24/29	14/23	2
SBS/55	34/35	34/35	14/23	2




DISTILLATION HEADS

Distillation head, splash, anti-foaming, sloping.


<i>Code No.</i>	<i>Cone to fit flask</i>	<i>Cone to fit condenser</i>	<i>Sockets size</i>	<i>Std. pack qty.</i>
SDH 10/23	24/29	19/26	14/23	1

Distillation head, splash, double bulb, sloping.




<i>Code No.</i>	<i>Cone to fit flask</i>	<i>Cone to fit condenser</i>	<i>Sockets size</i>	<i>Std. pack qty.</i>
SDH 14/23	24/29	19/26	14/23	1

Distillation heads, splash, sloping.



<i>Code No.</i>	<i>Cone to fit flask</i>	<i>Cone to fit condenser</i>	<i>Std. pack qty.</i>
SPS/22	19/26	19/26	2
SPS/23	24/29	19/26	2
SPS/33	24/29	24/29	1
SPS/44	29/32	29/32	1
SPS/55	34/35	34/35	1

Distillation heads, splash, vertical.



<i>Code No.</i>	<i>Cone to fit flask</i>	<i>Cone to fit condenser</i>	<i>Std. pack qty.</i>
SPV/22	19/26	19/26	5
SPV/23	24/29	19/26	5
SPV/33	24/29	24/29	2