


# EXTRACTORS



## Extractors, plain funnel distributors.

Code No.	Capacity ml	Socket size	Cone size	Std. pack qty.
Extractor body				
EJB 10/35	10	34/35	24/29	1
EJB 10/55	10	34/35	34/35	1



## Extractors, disc-baffle distribution.

Code No.	Capacity ml	Socket size	Cone size	Overall height mm	Std. pack qty.
Extractor					
ELH 6/35	60	34/35	24/29	250	1
Distributor					
EDHB/6					1



## Extractors, sintered disc distributors.

Code No.	Capacity ml	Socket size	Cone size	Overall height mm	Std. pack qty.
Extractor					
ELH 40/35	400	34/35	24/29	890	1
ELH 1/59	1000	55/44	34/35	650	1
Distributor					
EDHS/40	For ELH 40/35				1
EDHS/1	For ELH 1/59				1

## Extractors, liquid/liquid.


Code No.	Capacity ml	Socket size	Cone size	Overall height mm	Std. pack qty.
Extractor					
ELL 5/35	100	34/35	24/29	585	1
ELL 10/55	100	34/35	34/35	585	1
ELL 20/35	200	34/35	24/29	770	1



## Distributors, sintered.

Code No.	Type	To fit
EDLS 5	Sinter	ELL 5/35
EDLS 10	Sinter	ELL 10/35
EDLS 20	Sinter	ELL 20/35

## Distributors, disc baffle.

Code No.	
EDLD 5	For ELL 5/35
EDLD 10	For ELL 10/35
EDLO 20	For ELL 20/35





# EXTRACTORS

## Extractors, plain body.

Code No.	Capacity ml	Socket size	Cone size	Overall height mm	Std. pack qty.
<b>Body</b>					
EJB 10/35	60	34/35	24/29	200	1
EJB 10/55	60	34/35	34/35	200	1
EJB/37	100	45/40	24/29	230	1
EJB/35	190	34/35	24/29	350	2
EJB/39	200	55/44	24/29	250	1



## Extractors, Soxhlet.

Code No.	Nominal capacity (siphoning volume) ml	Socket or flange size	Cone size	Overall height mm	Std. pack qty.
ESC 6/21	6	19/26	14/23	150	1
ESC 30/32	20	24/29	19/26	175	1
ESC 40/43	40	29/32	24/29	200	1
ESC 60/53	60	34/35	24/29	205	2
ESC 60/55	60	34/35	34/35	205	2
ESC 10/55	100	34/35	34/35	240	2
ESC 10/63	100	40/38	24/29	220	2
ESC 20/55	200	34/35	34/35	360	1
ESC 20/83	200	50/42	24/29	250	1
ESC 60/5	600	75 mm	34/35	370	1