


JACKETED REACTION VESSEL

- * Buttress end inlet and outlet for oil or water circulation.
- * Special valve ensuring no "dead space".
- * Tangential entry and exit of circulating fluid.
- * 1 litre capacity.
- * Flat flange fitting for use with multi-neck lids.



Jacket Reaction Vessel.

Code No.	Capacity	Flange size	Buttress end size
JRV/1	1 litre	100 mm	16 mm

This apparatus was developed with the co-operation of Ciba Geigy (U.K.) Process Development Laboratories. We are pleased to acknowledge their invaluable assistance in its design.