


JOINTS


SCIENTIFIC-LAB. GLASS AB



Joints, cone, angle drip, double.

Code No.	Cone size	Shank length mm	O/D of shank mm	Min. bore mm	Wall mm	Std. pack qty.
CDD 65/B 14	14/23	65	14.0	6.5	1.2	5
CDD 190/B 19	19/26	190	18.0	8.5	1.8	5
CDD 150/B 24	24/29	150	23.0	14.5	2.0	5
CDD 120/B 29	29/32	120	28.0	19.5	2.0	5

Joints, cone, angle drip, single.



Code No.	Cone size	Shank length mm	O/D of shank mm	Min. bore mm	Wall mm	Std. pack qty.
CD 110/B 10	10/19	110	9.0	3.5	1.2	20
CD 110/B 14	14/23	110	14.0	6.5	1.2	20
CD 110/B 19	19/26	110	18.0	8.5	1.8	20
CD 110/B 24	24/29	110	23.0	14.5	2.0	20
CD 110/B 29	29/32	110	28.0	19.5	2.0	20
CD 110/B 34	34/35	110	33.5	22.5	2.0	20

Joints, cone, plain end, double.



Code No.	Cone size	Shank length mm	O/D of shank mm	Min. bore mm	Wall mm	Std. pack qty.
CND 50/B 10	10/19	50	9.0	5.0	1.3	10
CND 80/B 14	14/23	80	14.0	8.5	1.5	10
CND 80/B 19	19/26	80	18.0	12.0	1.8	10
CND 95/B 24	24/29	95	23.0	16.5	2.0	10
CND 100/B 29	29/32	100	28.0	21.0	2.0	10
CND 150/B 34	34/35	150	33.5	26.0	2.0	10

Joints, cone, plain end cut back at angle, double.



Code No.	Cone size	Shank length mm	O/D of shank mm	Min. bore mm	Wall mm	Std. pack qty.
CPD 60/B 10	10/13	60	9.0	6.0	1.2	5

Joints, cone, plain end, single.



Code No.	Cone size	Shank length mm	O/D of shank mm	Min. bore mm	Wall mm	Std. pack qty.
CN 110/B 7	7/16	110	7.0	2.5	1.2	20
CN 110/B 10	10/19	110	9.0	5.0	1.3	20
CN 110/B 12	12/21	110	12.0	6.5	1.5	20
CN 110/B 14	14/23	110	13.0	8.5	1.5	20
CN 110/B 19	19/26	110	18.0	12.0	1.8	20
CN 110/B 24	24/29	110	23.0	16.5	2.0	20
CN 110/B 29	29/32	110	28.0	21.0	2.0	20
CN 110/B 34	34/35	110	33.5	26.0	2.0	20
CN 110/B 40	40/38	110	39.0	30.5	2.4	20
CN 110/B 45	45/40	110	43.0	35.0	2.5	20
CN 110/B 50	50/42	110	49.0	39.0	2.8	20
CN 110/B 55	55/44	110	53.0	43.5	2.8	20


JOINTS

SCIENTIFIC-LAB.GLASS AB



Joins, cone, plain end, with hooks.

Code No.	Cone size	Shank length mm	O/D of shank mm	Min. bore mm	Wall mm	Std. pack qty.
CHB 110/B 10	10/19	110	9.0	5.0	1.3	5
CHB 110/B 14	14/23	110	14.0	8.5	1.5	5
CHB 110/B 19	19/26	110	18.0	12.0	1.8	5
CHB 110/B 24	24/29	110	23.0	16.5	2.0	5
CHB 110/B 29	29/32	110	28.0	21.0	2.0	5



Joins, cone with long tip, double.

Code No.	Cone size	Shank length mm	O/D of shank mm	Min. bore mm	Wall mm	Std. pack qty.
CBD 55/B 14	14/23	55	14.0	6.0	1.2	20
CBD 60/B 19	19/26	60	18.0	8.0	1.8	20
CBD 95/B 24	24/29	95	23.0	14.0	2.0	10
CBD 125/B 29	29/32	125	28.0	19.0	2.0	10
CBD 165/B 34	34/35	165	33.5	22.0	2.0	10



Joins, cone with long tip, single.

Code No.	Cone size	Shank length mm	O/D of shank mm	Min. bore mm	Wall mm	Std. pack qty.
CB 90/B7	7/16	90	7.0	2.0	1.2	20
CB 90/B 10	10/19	90	9.0	3.0	1.5	20
CB 90/B 14	14/23	90	14.0	6.0	1.2	20
CB 90/B 19	19/26	90	18.0	8.0	1.8	20
CB 90/B 24	24/29	90	23.0	14.0	2.0	20
CB 90/B 29	29/32	90	28.0	19.0	2.0	20
CB 90/B 34	34/35	90	33.5	22.0	2.0	20
CB 90/B 40	40/38	90	39.0	27.0	2.2	10
CB 90/B 45	45/40	90	43.0	31.0	2.5	10
CB 90/B 55	55/44	90	53.0	40.0	2.8	10



Joins, cone with stem, reduced shank, single.

Code No.	Cone size	Shank length mm	O/D Of shank mm	Wall mm	Stem length mm	Min. bore mm	Std. pack qty.
CS 45/B 10	10/19	45	9.0	1.2	95	2.5	10
CS 45/B 14	14/23	45	7.0	1.2	150	3.5	10
CS 60/B 19	19/26	60	10.0	1.5	150	6.5	10
CS 70/B 24	24/29	70	10.0	1.5	150	6.5	10
CS 70/B 29	29/32	70	10.0	1.5	150	6.4	1





JOINTS

Joins, flat flange (shank length 110 mm).

Code No.	Bore mm	O/D of shank mm	Wall mm	O/D of flange mm	Std. pack qty.
G 1/2	2	8	2.2	20	1
G 1/5	5	10	2.2	25	1
G 1/10	10	15	2.2	30	1
G 1/15	15	21	3.5	45	1
G 1/25	25	32	2.2	55	1
G 1/35	35	40	2.7	75	1
G 1/50	50	55	3.5	85	1
G 1/75	75	85	5.0	120	1
G 1/100	100	110	5.0	150	1

Joins, screwthread



Code No.	Thread size	Approx. shank length mm	Shank O/D mm	Min. bore mm	Std. pack qty.
JS 29/3	13	110	11.5	7.5	20
JS 29/4	18	110	15.5	10.5	20
JS 29/5	24	110	22.5	15.5	20
JS 29/6	28	110	24.5	19.7	20

FOR FITTINGS SEE PAGE 37

Joins, spherical, ball and cup (shank length 100 mm, except FS 51 which is 120 mm)



Code No.	Joint size	O/D Of shank mm	Wall mm	Bore mm	Std. pack qty.
Ball joint (male)					
J 14/S 13C	S 13	9	2.0	2	10
J 14/S 13	S 13	9	1.7	5	10
J 14/S 19	S 19	14	1.7	9	10
J 14/S 29	S 29	22	2.2	15	10
J 14/S 35	S 35	30	2.2	21	10
J 14/S 41	S 41	34	2.2	27	5
J 14/S51	S 51	43	2.6	35	5
Cup joint (female)					
J 15/S 13C	S 13	9	2.0	2	10
J 15/S 13	S 13	9	1.7	5	10
J 15/S 19	S 19	14	1.7	9	10
J 15/S 29	S 29	22	2.2	15	10
J 15/S 35	S 35	30	2.2	21	5
J 15/S 41	S 41	34	2.2	27	5
J 15/S 51	S 51	43	2.6	35	5


Joins, socket, bulged shoulder, single (shank length 110 mm).



Code No.	Socket size	O/D of shank mm	Min. bore mm	Wall mm	Std. pack qty.
SB 110/B 7	7/16	7.0	3.0	1.0	2
SB 110/B 10	10/19	10.0	6.0	1.0	5
SB 110/B 12	12/21	12.0	7.5	1.5	2
SB 110/B 14	14/23	14.0	9.5	1.5	5
SB 110/B 19	19/26	19.0	13.0	1.8	2
SB 110/B 24	24/29	23.0	17.5	2.0	2
SB 110/B 29	29/32	28.0	22.5	2.0	5
SB 110/B 34	34/35	33.5	27.5	2.0	2

JOINTS

SCIENTIFIC-LAB. GLASS AB



Joins, socket, rolled shoulder, double.

Code No.	Socket size	Shank length mm	O/D of shank mm	Min. bore mm	Wall mm	Std. pack qty.
SRD 85/B 7	7/16	85	11.0	3.5	1.5	2
SRD 85/B 10	10/19	85	14.0	8.0	1.5	20
SRD 85/B 14	14/23	85	18.0	12.0	1.8	20
SRD 85/B 19	19/26	85	23.0	16.0	2.0	20
SRD 145/B 24	24/29	145	28.0	21.0	3.0	5
SRD 75/B 29	29/32	75	33.5	26.0	2.0	10
SRD 100/B 34	34/35	100	40.0	31.0	2.2	2

Joins, socket, rolled shoulder, single.




Code No.	Socket size	Shank length mm	O/D of shank mm	Min. bore mm	Wall mm	Std. pack qty.
SR 110/B 7	7/16	110	11.0	5.5	1.5	20
SR 110/B10	10/19	110	14.0	8.0	1.5	20
SR 110/B12	12/21	110	15.5	10.0	1.5	20
SR 110/B 14	14/12	110	18.0	12.0	1.8	20
SR 110/B 19	19/26	110	23.0	16.0	2.0	20
SR 110/B 24	24/29	110	28.0	21.0	2.0	20
SR 110/B 29	29/32	110	33.5	26.0	2.0	20
SR 110/B 34	34/35	110	40.0	31.0	2.2	20
SR 175/B 34	34/35	175	40.0	31.0	2.2	10
SR 110/B 40	40/38	110	46.0	36.0	2.4	20
SR 110/B 45	45/40	110	53.0	41.0	2.8	20
SR 110/B 50	50/42	110	57.0	45.5	2.8	10
SR 110/B 55	55/44	110	64.0	50.0	3.0	10

Joins, socket, rolled shoulder, single with hooks (shank length 110 mm).



Code No.	Socket size	O/D of shank mm	Min. bore mm	Wall mm	Std. pack qty.
SH 110/B 14	14/23	18.0	12.0	1.8	2
SH 110/B 19	19/26	23.0	16.0	2.0	2
SH 110/B 24	24/29	28.0	21.0	2.0	2
SH 110/B 29	29/32	33.5	26.0	2.0	2

Joins, socket with cup (shank length 75 mm).



Code No.	Socket size	O/D of shank mm	Shank bore mm	Wall mm	Cup bore mm	Std. pack qty.
SC 75/B 10	10/19	14.0	10	1.5	15	5
SC 75/B 14	14/23	18.0	13	1.8	20	5
SC 75/B 19	19/26	23.0	18	2.0	25	5
SC 75/B 24	24/29	28.0	23	2.0	30	5
SC 75/B 29	29/32	33.5	26	2.0	35	5