



SCREW THREADS & FITTINGS

Pipette weighing, Lunge-Rey

Code No.		Joint size	Capacity ml
LR 1	Pipette	10/19	4
TT 4/0	Cover tube	10/19	4
FC 25/B10	Cover flask	10/19	25



Receiver with vent.

Code No.	Nominal capacity ml	Socket size	Approx. overall height mm	Std. pack qty.
RV 1	5	10/19	95	10
RV 2	10	14/23	130	10
RV 3	25	14/23	70	5

Screwthread Joint Fittings (Plastic caps, silicone rubber rings and PTFE washers)



Thread size	Acceptable dia. through screwcap mm	Plastics cap code no.	Rubber ring code no.	PTFE washer code no.
13	6.0— 7.0	JC 13/3	JR 13/3	JW 13/3
18	7.0— 8.5	JC 18/3	JR 18/3	JW 18/3
24	6.0— 7.0	JC 24/3	JR 24/3	JW 24/3
24	8.5— 9.5	JC 24/4	JR 24/4	JW 24/4
24	10.5—11.5	JC 24/5	JR 24/5	JW 24/5
28	6.0— 7.0	JC 28/3	JR 28/3	JW 28/3
28	10.5—11.5	JC 28/4	JR 28/4	JW 28/4
28	18.0—19.0	JC 28/5	JR 28/5	JW 28/5

Plastics cap



Code No.	Std. pack qty.
JC 13/3	50
JC 18/3	50
JC 24/3	50
JC 28/3	50
JC 28/4	50
JC 28/5	50

PTFE washer

Code No.	Std. pack qty.
JW 13/3	100
JW 18/3	100
JW 24/3	10
JW 24/4	10
JW 24/5	20
JW 28/3	10
JW 28/4	10
JW 28/5	10



Rubber ring

Code No.	Std. pack qty.
JR 13/3	100
JR 18/3	50
JR 24/3	20
JR 24/4	10
JR 24/5	20
JR 28/3	5
JR 28/4	20
JR 28/5	20

