


# TUBES


## Tubes, test.

Code No.	Socket size	Approx. overall length mm	Max. body diam. mm	Approx. capacity ml	Std pack qty.
TT6/0	10/19	65	14	2	20
TT10/0	10/19	100	14	5	20
TT10/12	12/21	100	17	10	10
TT10/1	14/23	100	18	11	20
TT12/1	14/23	125	18	15	20
TT15/1	14/23	150	18	19	20
TT10/2	19/26	100	23	19	20
TT15/2	19/26	150	23	31	20
TT15/3	24/29	150	29	54	20
TT20/3	24/29	200	29	75	10
TT25/5	34/35	250	41	195	5



## Tubes, 'U', with cones.

Code No.	Cone size	Std. pack qty.
TU/1	14/23	1
TU/2	19/28	1



## Tubes, 'U', with sockets.

Code No.	Socket size	Std. pack qty.
TUS/1	14/23	2
TUS/2	19/26	2



## Non-return valves, plain.


NVP/1

Other designs of non-return valve, to customers own requirements, are made to order.



## Tubes, drying, angled.

Code No.	Cone size	Std. pack qty.
DTS/0	10/19	5
DTS/1	14/23	5
DTS/2	19/26	5
DTS/3	24/29	5
DTS/4	29/32	5



# TUBES

## Tubes, bent, with cone.

<i>Code No.</i>	<i>Cone size</i>	<i>Stem length mm</i>	<i>Std. pack qty.</i>
MCS/0B	10/19	100	2
MCS/1B	14/23	150	5
MCS/2B	19/26	150	5
MCS/3B	24/29	150	5


49



## Tubes/air inlets, straight with cone.

<i>Code No.</i>	<i>Stem length mm</i>	<i>Cone Size</i>	<i>Std. pack qty.</i>
MCS 0	100	10/19	5
MCS 1	150	14/23	5
MCS 2	150	19/26	5
MCS 3	150	24/29	5
MCS 4	150	24/29	5

# THERMOMETER POCKETS



## Thermometer pockets.

<i>Code No.</i>	<i>Cone size</i>	<i>Stem length mm</i>	<i>Std. pack qty.</i>
TPS 1	14/23	45	5
TPS 2	19/26	200	5